

BLACK NYLON SOFT BAG SPENGLER, LARGE

Ref. 535019



Description

Nylon pouch 25x20x5 cm with éclair® zip closure

 07/04/2026 02:20:56 – For information purposes only, can be subject to change

BLACK NYLON SOFT BAG SPENGLER, LARGE

Ref. 535019



Features

CE Marking	NONE	Medical device	No
Bar code	3700446002008	HS Code	90330000
Warranty (months)	60	Country	China
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.111
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

📅 07/04/2026 02:20:56 – For information purposes only, can be subject to change