

LARGE ADULT 31-40CM X20

Ref. DSU231067

SUNTECH



CE



Single-use large adult cuffs for comfortable and accurate blood pressure monitoring.

Description

Set of 20 Suntech "SOFT" single-patient ambands, burgundy.
Adult wide size 31-40cm, rectus bayonet connection.

06/04/2026 23:39:36 – For information purposes only, can be subject to change

LARGE ADULT 31-40CM X20

Ref. DSU231067

SUNTECH

Features

CE Marking	CE	Medical device	Class I
Bar code	3700446014186	HS Code	9018905040
Country	China		
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	24	Packed unit / Height (cm)	24
Packed unit / Length (cm)	24	Packed unit / Weight	2.004
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

📅 06/04/2026 23:39:36 – For information purposes only, can be subject to change