

TIGHT 38-50CM X20

Ref. DSU231069

SUNTECH



CE



07/04/2026 02:14:51 – For information purposes only, can be subject to change

TIGHT 38-50CM X20

Ref. DSU231069

SUNTECH

Extra-large single-use cuffs – adult size (23–33 cm), Bayonet connector, designed for specific clinical needs.

Description

Set of 20 Suntech "SOFT" single-patient armbands, brown.
 Bib shorts size 38-50cm, rectus bayonnette connection.

Features

CE Marking	CE	Medical device	Class I
Bar code	3700446014407	HS Code	9018905040
Country	China		

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	2.49

SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

OTHER FEATURES

Packed unit / Quantity	1		
------------------------	---	--	--

📅 07/04/2026 02:14:51 – For information purposes only, can be subject to change