

**SET OF 5 SUNTECH MONOBLOC REUSABLE
ARMBANDS ADULT PLUS (28-40CM) -
BAYONET METAL CONNECTOR**

SUNTECH

Ref. DSU231114



CE



07/04/2026 01:18:11 – For information purposes only, can be subject to change

SET OF 5 SUNTECH MONOBLOC REUSABLE ARMBANDS ADULT PLUS (28-40CM) - BAYONET METAL CONNECTOR

SUNTECH

Ref. DSU231114

Set of 5 Suntech Adulte Plus cuffs (28–40 cm), reusable with metal Bayonet connector for reliable and long-lasting measurement.

Description

Reusable nylon cuff sold in boxes of 5 units.
Single tubing, bayonette connection 28-40cm.

Features

CE Marking	CE	Medical device	Class I
Bar code	3661041008829	HS Code	9018905040
Warranty (months)	24	Country	China

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	25	Packed unit / Height (cm)	10
Packed unit / Length (cm)	23	Packed unit / Weight	0.6

SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

OTHER FEATURES

Packed unit / Quantity	1
------------------------	---

07/04/2026 01:18:11 – For information purposes only, can be subject to change