

**MIXTER-O'SHAUGNESSY ARTERY FORCEPS**  
**90° 16 CM**

Ref. IP49016



CE

DISPOSITIF MEDICAL  
CLASSE  
**I**  
MEDICAL DEVICE

GARANTIE  
2 ANS  
**2**  
2 YEARS  
GUARANTEE

STAINLESS  
**INOX**  
STEEL

07/04/2026 02:36:07 – For information purposes only, can be subject to change

**MIXTER-O'SHAUGNESSY ARTERY FORCEPS**  
**90° 16 CM**

Ref. IP49016



**Features**

CE Marking	CE	Medical device	Class I
Bar code	3661040138305	HS Code	90189084
Warranty (months)	24	Country	Pakistan
<b>DIMENSIONS AND WEIGHT</b>			
Packed unit / Width (cm)	9	Packed unit / Height (cm)	1
Packed unit / Length (cm)	22	Packed unit / Weight	0.07
Bare Unit / Weight (kg)	0.07		
<b>SPECIFIC STORAGE CONDITIONS</b>			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
<b>OTHER FEATURES</b>			
Packed unit / Quantity	1		

📅 07/04/2026 02:36:07 – For information purposes only, can be subject to change