

**MIXTER-O'SHAUGNESSY ARTERY FORCEPS**  
**90° 18 CM**

**Ref. IP49018**



**Features**

CE Marking	CE	Medical device	Class I
Bar code	3661040138312	HS Code	90189084
Warranty (months)	24	Country	Pakistan
<b>DIMENSIONS AND WEIGHT</b>			
Packed unit / Width (cm)	8	Packed unit / Height (cm)	2
Packed unit / Length (cm)	19	Packed unit / Weight	0.07
Bare Unit / Weight (kg)	0.07		
<b>SPECIFIC STORAGE CONDITIONS</b>			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
<b>OTHER FEATURES</b>			
Packed unit / Quantity	1		

06/04/2026 23:32:13 – For information purposes only, can be subject to change