


BAG EIFFEL, DE BOISSY

Ref. MDEIFF01

De Boissy
— DEPUIS 1982



📅 07/04/2026 03:34:22 – For information purposes only, can be subject to change

BAG EIFFEL, DE BOISSY

Ref. MDEIFF01



DEPUIS 1982

The chic and simple large volume bag, easy to organize

Description

Chic and simple at the same time, the large EIFFEL bag recalls the shape of the traditional leather bag of doctors. Thanks to its very large bidirectional opening and multiple pockets/compartments, it makes it easy to store and reach all the care equipment and diagnostic instruments that health professionals need.

Of a mottled gray, the polyester fabric is resistant, durable, water-repellent and therefore washable, it is accompanied by a reflective strip for safety, and YKK® slides for their quality.

- The large volume interior main space includes 8 pockets, 2 removable velcro dividers, 1 zipped mesh pouch. With a very large bidirectional zip opening, it gives very easy access to diagnostic instruments and care equipment
- A ventral pocket (40x20 cm) opens completely with on one side 2 zipped storage nets and 2 flat pockets, and on the other side several elastic loops that allow to store small care supplies: pliers, scissors, pen lamp, thermometer, bandages, pens, etc...
- Four external zipped pockets, one for A4 documents and another on the side with elastic to slide the bottle of hydroalcoholic gel which remains easily accessible at any time
- Padded and solid handles, a padded and adjustable shoulder strap - A name tag to identify the bag
- A non-slip bottom to reduce friction and wear
- Dimensions w x h x d: 40 x 35 x 24 cm- Weight: 1.3 kg
- Maximum recommended content: 6 kg / 24 liters
- Sold empty

📅 07/04/2026 03:34:22 – For information purposes only, can be subject to change

BAG EIFFEL, DE BOISSY

Ref. MDEIFF01



Features

CE Marking	NONE	Medical device	No
Bar code	3661040218939	HS Code	420229000
Warranty (months)	12	Country	China

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	37	Bare Unit / Width (cm)	35
Packed unit / Height (cm)	10	Bare Unit / Height (cm)	24
Packed unit / Length (cm)	45	Bare Unit / Length (cm)	40
Packed unit / Weight	nc	Bare Unit / Weight (kg)	1.29

SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

OTHER FEATURES

Packed unit / Quantity	1
------------------------	---

07/04/2026 03:34:22 – For information purposes only, can be subject to change