


## MINI BAG EIFFEL, DE BOISSY

Ref. MDEIFF02

*De Boissy*  
— DEPUIS 1982



08/04/2026 01:32:08 – For information purposes only, can be subject to change

## MINI BAG EIFFEL, DE BOISSY

Ref. MDEIFF02



DEPUIS 1982

### The mini version of the EIFFEL Medical Bag

## Description

Chic and simple at the same time, the EIFFEL mini bag recalls the shape of the traditional leather bag of doctors. Thanks to its bidirectional opening and multiple pockets/compartments, it makes it easy to store and reach all the care equipment and diagnostic instruments that healthcare professionals need. Of a mottled gray, the polyester fabric is resistant, durable, water-repellent and therefore washable, it is accompanied by a reflective strip for safety, and YKK® slides for their quality.

- A main space that includes 6 pockets and 1 removable velcro divider. The large bidirectional zip opening gives very easy access to diagnostic instruments and care equipment
- A ventral pocket (31x15cm) opens completely with on one side 2 zipped storage nets and 2 flat pockets, and on the other side several elastic loops that allow to store small care supplies: pliers, scissors, pen lamp, thermometer, bandages, pens, etc...
- There are also four small external zipped pockets and another on the side with elastic to slide the bottle of hydroalcoholic gel that remains easily accessible at any time
- Padded and solid handles, a padded and adjustable shoulder strap - A name tag to identify the bag
- A non-slip bottom to reduce friction and wear
- Dimensions w x h x d: 32 x 28 x 20 cm- Weight: 0.95 kg
- Maximum recommended content: 4 kg / 15 liters
- Sold empty

📅 08/04/2026 01:32:08 – For information purposes only, can be subject to change

# MINI BAG EIFFEL, DE BOISSY

Ref. MDEIFF02



## Features

CE Marking	NONE	Medical device	No
Bar code	3661040218946	HS Code	4202290000
Warranty (months)	12	Country	China

### DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	28	Bare Unit / Width (cm)	28
Packed unit / Height (cm)	8	Bare Unit / Height (cm)	20
Packed unit / Length (cm)	42	Bare Unit / Length (cm)	32
Packed unit / Weight	nc	Bare Unit / Weight (kg)	0.94

### SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

### OTHER FEATURES

Packed unit / Quantity	1
------------------------	---

08/04/2026 01:32:08 – For information purposes only, can be subject to change