

FOOT STOOL DOUBLE STEP STAINLESS STEEL

Ref. OS0050



07/04/2026 03:24:06 – For information purposes only, can be subject to change

FOOT STOOL DOUBLE STEP STAINLESS STEEL

Ref. OS0050



Ideal for medical practices, infirmaries, treatment rooms...

Description

- Stainless steel steps
- Steps covered with non-slip plastic coating
- Maximum load: 150 kg

Features

CE Marking	CE	Medical device	Class I
Bar code	3661040162614	HS Code	9402900000
Warranty (months)	12	Country	China

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	35	Bare Unit / Width (cm)	24.5
Packed unit / Height (cm)	12	Bare Unit / Height (cm)	34
Packed unit / Length (cm)	57	Bare Unit / Length (cm)	40
Packed unit / Weight	5.7	Bare Unit / Weight (kg)	5.1

SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

OTHER FEATURES

Packed unit / Quantity	1
------------------------	---

📅 07/04/2026 03:24:06 – For information purposes only, can be subject to change