

STAINLESS STEEL EXAMINATION COUCH (QUICK ASSEMBLING) BLACK

Ref. OSA003_N



06/04/2026 23:11:13 – For information purposes only, can be subject to change

STAINLESS STEEL EXAMINATION COUCH (QUICK ASSEMBLING) BLACK

Ref. OSA003_N



Quick Mount System

Description

- Rapid assembly system
- Square stainless steel tube frame
- Adjustable headboard with automatic high resistance rack
- Maximum weight: 150 kg
- Headboard and seat covered with 5 cm thick very high density foam (30 kg/m³), covered by fire-resistant washable plasticised fabric (fire resistance class: M1)
- Removable roll holder (for 50 cm wide rolls) that can be located at head or feet

OPTIONS

Réf. OTA0047 Heel support without clamps (pair)

Réf. OTA0041 Leg holders without clamps (pair)

Réf. OTA0044 Fixing clamps (pair)

 06/04/2026 23:11:13 – For information purposes only, can be subject to change

STAINLESS STEEL EXAMINATION COUCH (QUICK ASSEMBLING) BLACK



Ref. OSA003_N

Features

CE Marking	CE	Medical device	Class I
Bar code	3661040165844	HS Code	9402900000
Warranty (months)	36	Country	China
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	15	Bare Unit / Width (cm)	63
Packed unit / Height (cm)	68	Bare Unit / Height (cm)	70
Packed unit / Length (cm)	190	Bare Unit / Length (cm)	190
Packed unit / Weight	27.5	Bare Unit / Weight (kg)	24.5
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

06/04/2026 23:11:13 – For information purposes only, can be subject to change