

**INSTRUMENT TROLLEY S.S. 2 SHELVES
WITH SIDE RAILS 70X50X80 CM**

Ref. OSA7050



07/04/2026 00:12:47 – For information purposes only, can be subject to change

INSTRUMENT TROLLEY S.S. 2 SHELVES WITH SIDE RAILS 70X50X80 CM



Ref. OSA7050

High quality pedestal table, perfect for storing the necessary surgical instruments or consumables

Description

- Round tube structure in non-magnetic stainless steel for a better grip
- Stainless steel trays
- Removable galleries- Silent rubber wheels
- Permissible load: 80 kg
- Wheel diameter: 75 mm with brakes

Features

CE Marking	CE	Medical device	Class I
Bar code	3661040168319	HS Code	9402900000
Warranty (months)	12	Country	China

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	53	Bare Unit / Width (cm)	50
Packed unit / Height (cm)	11	Bare Unit / Height (cm)	80
Packed unit / Length (cm)	85	Bare Unit / Length (cm)	70
Packed unit / Weight	11.5	Bare Unit / Weight (kg)	11.5

SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

OTHER FEATURES

Packed unit / Quantity	1
------------------------	---

 07/04/2026 00:12:47 – For information purposes only, can be subject to change