

KIT EXAMINATION COUCH + STEPPER 1 STEP- EPOXY-BLACK - HOLTEX

Ref. OSKIT003



07/04/2026 02:14:32 – For information purposes only, can be subject to change

KIT EXAMINATION COUCH + STEPPER 1 STEP- EPOXY-BLACK - HOLTEx



Ref. OSKIT003

Full medical office

Description

Delivered in a single carton

Set composed of:- 1

sofa black color (ref OSA003E_N)- 1 step

1 step (OS0051E)

The sofa with structure in square tube in white epoxy lacquered steel:

-Adjustable headrest 75 ° with 2 self-locking mechanical racks.

-Head and seat padded with very high density foam (30 kg / m³), whose thickness is 5 cm, covered with washable plasticized fabric, fireproof (fire resistance classification: M1)

-Removable roller holder (for rolls of 50 cm wide) positionable at the head or feet

-Maximum load: 150 kg

The step 1 steps:

-white epoxy lacquered steel and non-slip plastic coating

- Load: 150 kg

Stainless steel or Epoxy: what's the difference?

Epoxy: the economic choice. The steel is simply coated with a resin to delay corrosion.

Stainless steel: it is steel modified to become completely STAINLESS steel, it does not rust.

And did you know?

Today, more than 90% of stainless steel is made from recycled steel, it's good for the environment!

 07/04/2026 02:14:32 – For information purposes only, can be subject to change

KIT EXAMINATION COUCH + STEPPER 1 STEP- EPOXY-BLACK - HOLTEX



Ref. OSKIT003

Features

CE Marking	CE	Medical device	Class I
Bar code	3661040168395	HS Code	9402900000
Warranty (months)	24	Country	China
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	15	Packed unit / Height (cm)	68
Packed unit / Length (cm)	190	Packed unit / Weight	34
Bare Unit / Weight (kg)	31.5		
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

📅 07/04/2026 02:14:32 – For information purposes only, can be subject to change