

# ALUMINIUM FOLDABLE STRETCHER IN 2 (LENGHT AND LARG) + POUCH, ORANGE

Ref. SBRAN08\_O



07/04/2026 00:44:39 – For information purposes only, can be subject to change

## ALUMINIUM FOLDABLE STRETCHER IN 2 (LENGHT AND LARG) + POUCH, ORANGE

Ref. SBRAN08\_O



**Folding stretcher with balance wheel system**

### Description

Folding stretcher with balancing system

- High resistance aluminium structure and waterproof fabric, fire resistant and anti-cracking
- This balancing stretcher folds in length and width
- Lightweight, easy to carry
- Dimensions (L x W x H): 208 X 55 x 13 cm, unfolded;  
104 x 17 x 9 cm, folded
- Weight: 8 kg
- Maximum load carried: 159 kg
- Single colour: orange
- Delivered with a stretcher bag

 07/04/2026 00:44:39 – For information purposes only, can be subject to change

# ALUMINIUM FOLDABLE STRETCHER IN 2 (LENGHT AND LARG) + POUCH, ORANGE



Ref. [SBRAN08\\_O](#)

## Features

CE Marking	CE	Medical device	Class I
Bar code	3661040094632	HS Code	94029000
Warranty (months)	12	Country	China
<b>DIMENSIONS AND WEIGHT</b>			
Packed unit / Width (cm)	20	Bare Unit / Width (cm)	54.5
Packed unit / Height (cm)	10	Bare Unit / Height (cm)	13.5
Packed unit / Length (cm)	10	Bare Unit / Length (cm)	208
Packed unit / Weight	5.6	Bare Unit / Weight (kg)	4.7
<b>SPECIFIC STORAGE CONDITIONS</b>			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
<b>OTHER FEATURES</b>			
Packed unit / Quantity	1		

07/04/2026 00:44:39 – For information purposes only, can be subject to change