

PE/PLASTIC SCOOP STRETCHER, YELLOW

Ref. SBRAN13_J



07/04/2026 00:48:38 – For information purposes only, can be subject to change

PE/PLASTIC SCOOP STRETCHER, YELLOW

Ref. SBRAN13_J



Description

- «Spoon» function performed by the two longitudinal blades, to position the stretcher under the patient without moving him. The scoop stretcher can be used to lift, move, or transport patients.
- Length adjustable on 4 levels according to the patient's size, from 162 to 234 cm
- High resistance aluminium structure
- Dimensions (L x W x H): 162-234 x 44 x 7 cm
- Weight: 7 kg
- Maximum load carried: 220 kg

PE/PLASTIC SCOOP STRETCHER, YELLOW

Ref. SBRAN13_J



Features

CE Marking	CE	Medical device	Class I
Bar code	3661040094670	HS Code	94029000
Warranty (months)	12	Country	China
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	47	Bare Unit / Width (cm)	45
Packed unit / Height (cm)	8	Bare Unit / Height (cm)	7
Packed unit / Length (cm)	165	Bare Unit / Length (cm)	162
Packed unit / Weight	10.8	Bare Unit / Weight (kg)	8
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

07/04/2026 00:48:38 – For information purposes only, can be subject to change