

ALUMINIUM ALLOY STAIR STRETCHER - W 38 CM - 2 WHEELS, ORANGE

Ref. SBRAN14_O



06/04/2026 22:41:14 – For information purposes only, can be subject to change

ALUMINIUM ALLOY STAIR STRETCHER - W 38 CM - 2 WHEELS, ORANGE

Ref. SBRAN14_O



Ideal for transferring a patient, going up and down stairs

Description

- High strength aluminium alloy structure
- Ideal for transferring a patient, climbing up and down the stairs
- Two fixed rear wheels allow movement on most surfaces with minimal effort
- The safety of patients during the transport process is ensured by 3 restraining straps with quick release buckle
- At the front of the stretcher, two retractable arms
- Dimensions (L x W x H): 62 X 51 x 91 cm, unfolded;
90 x 17 x 59 cm, folded
- Weight: 8 kg
- Maximum load carried: 159 kg
- Soft seat
- Seat height: 49 cm
- Single colour: orange

 06/04/2026 22:41:14 – For information purposes only, can be subject to change

ALUMINIUM ALLOY STAIR STRETCHER - W 38 CM - 2 WHEELS, ORANGE



Ref. SBRAN14_O

Features

CE Marking	CE	Medical device	Class I
Bar code	3661040094687	HS Code	94029000
Warranty (months)	12	Country	China

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	56	Bare Unit / Width (cm)	52
Packed unit / Height (cm)	16	Bare Unit / Height (cm)	92
Packed unit / Length (cm)	98	Bare Unit / Length (cm)	84
Packed unit / Weight	10	Bare Unit / Weight (kg)	9

SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

OTHER FEATURES

Packed unit / Quantity	1
------------------------	---

06/04/2026 22:41:14 – For information purposes only, can be subject to change