

ALUMINIUM ALLOY STAIR STRETCHER - W 46 CM - 4 WHEELS, ORANGE

Ref. SBRAN15_O



07/04/2026 03:57:24 – For information purposes only, can be subject to change

ALUMINIUM ALLOY STAIR STRETCHER - W 46 CM - 4 WHEELS, ORANGE

Ref. SBRAN15_O



Ideal for transferring a patient, going up and down stairs with minimal effort

Description

- High strength aluminium alloy structure
- Equipped with armrests
- Ideal for transferring a patient, climbing up and down the stairs
- Back wheels are fixed and have brakes, the front two are swivel wheels, no brake. This allowing movement on most surfaces with a minimum of effort.
- The safety of patients during the transport process is ensured by 2 restraining straps with quick release buckle
- At the front of the stretcher, two retractable arms
- Unfolded dimensions: H98x D66.5
- Handle: 65.5 to 97 cm
- Dimensions folded : H98 x D60 x W20cm
- Rigid seat, 51 cm (deep) by 43 cm (wide)
- Seat height: 53 cm
- Weight: 10.5 kg
- Maximum load carried: 159 kg
- Rigid seat
- Single colour: orange

 07/04/2026 03:57:24 – For information purposes only, can be subject to change

ALUMINIUM ALLOY STAIR STRETCHER - W 46 CM - 4 WHEELS, ORANGE



Ref. SBRAN15_O

Features

CE Marking	CE	Medical device	Class I
Bar code	3661040094694	HS Code	94029000
Warranty (months)	12	Country	China

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	64	Bare Unit / Width (cm)	60
Packed unit / Height (cm)	20	Bare Unit / Height (cm)	98
Packed unit / Length (cm)	10	Bare Unit / Length (cm)	84
Packed unit / Weight	16	Bare Unit / Weight (kg)	14

SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

OTHER FEATURES

Packed unit / Quantity	1
------------------------	---

07/04/2026 03:57:24 – For information purposes only, can be subject to change