


CONVERTIBLE BACK PACK EIFFEL, DE BOISSY

Ref. MDEIFF03

De Boissy
— DEPUIS 1982



07/04/2026 00:54:43 – For information purposes only, can be subject to change

CONVERTIBLE BACK PACK EIFFEL, DE BOISSY

Ref. MDEIFF03

De Boissy
— DEPUIS 1982

It is worn in a backpack, hand or shoulder; A real Swiss army knife!

Description

Of a mottled gray, the polyester fabric is resistant, durable, water-repellent and therefore washable, it is accompanied by a reflective strip for safety, and YKK® slides for their quality. - Its main space

- is modular and contains 4 zipped compartments, extractable (fixed with velcro): 3 transparencies to quickly identify the contents, and 1 isothermal, supplied with its reusable cold gel pocket, to slip up to 20 bulbs
- An extension zip makes it possible to enlarge this main compartment to store more equipment
 - A second padded space can store computer / tablet and documents
 - This model also has the large ventral pocket (45x21 cm) which opens completely with on one side 2 zipped storage nets and 2 flat pockets, and on the other side several elastic loops that allow to store small care supplies: pliers, scissors, pen lamp, thermometer, bandages, pens, etc ...
 - Two external zipped pockets, plus one with elastic dedicated to the bottle of hydroalcoholic gel
 - The shoulder straps, handles and shoulder strap are padded, solid and concealable to adopt the wear that suits best
 - Dimensions w x h x d: 30 x 45 x 15/19 cm
 - Weight: 2.1 kg
 - Maximum recommended content: 6 kg / 19-25 liters
 - Sold empty

 07/04/2026 00:54:43 – For information purposes only, can be subject to change

CONVERTIBLE BACK PACK EIFFEL, DE BOISSY

Ref. MDEIFF03



Features

CE Marking	NONE	Medical device	No
Bar code	3661040218953	HS Code	4202129900
Warranty (months)	12	Country	China
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	50	Bare Unit / Width (cm)	15
Packed unit / Height (cm)	35	Bare Unit / Height (cm)	45
Packed unit / Length (cm)	50	Bare Unit / Length (cm)	30
Packed unit / Weight	2.1	Bare Unit / Weight (kg)	2.1
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

📅 07/04/2026 00:54:43 – For information purposes only, can be subject to change