

VITALIMATE MINI PATIENT MONITOR SPO2 SPENGLER

Ref. MON001



CE
0123

DISPOSITIF MEDICAL
CLASSE
II_b
MEDICAL DEVICE

GARANTIE
1 YEAR
GUARANTEE



07/04/2026 01:28:37 – For information purposes only, can be subject to change

VITALIMATE MINI PATIENT MONITOR SPO2 SPENGLER



Ref. MON001

Vitalimate MINI the polyvalent portable monitor for all health care settings

Description

Vitalimate MINI measures SpO2 and pulse rate. It is suitable for all patients types: adult, pediatric and neonat.

- High definition 3,5" (420x320 pixels). Can display figures and waveform
- A double function base. It allows up to 12 hours in continuous use or up to 24 hours if plugged on the base. The base also serves to store the device.
- Audible and visual configurable alarms
- Compact and easy to carry around: 280g, dimensions 175 x 85 x 32 mm

Box content: a monitor, a base, a rigid SpO2 clip adult sensor, a mains cord and a manual

Other sizes of SpO2 sensors are available and sold separately:

- Adult reusable clip SpO2 sensor : reference MONI010
- Pediatric reusable clip SpO2 sensor : reference MONI010B
- Neonatal reusable clip SpO2 sensor : reference MONI010D
- Reusable ear clip SpO2 sensor : reference MONI010E
- Adult silicon SpO2 sensor : reference MONI010A
- Pediatric silicon SpO2 sensor : reference MONI010C

VITALIMATE MINI PATIENT MONITOR SPO2

SPENGLER



Ref. MON001

Features

CE Marking	CE0123	Medical device	Class IIb
Bar code	6948309002025	HS Code	9018193010
Warranty (months)	12	Country	China
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	9	Packed unit / Height (cm)	33
Packed unit / Length (cm)	22	Packed unit / Weight	1.33
Bare Unit / Weight (kg)	0.585		
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	-20	Max temperature	55
Minimum humidity	10	Maximum humidity	93
OTHER FEATURES			

📅 07/04/2026 01:28:37 – For information purposes only, can be subject to change