

## EAR REUSABLE CLIP SPO2 SENSOR FOR VITALIMATE MINI

Ref. MONI010E



CE  
0123

DISPOSITIF MEDICAL  
CLASSE  
**II<sub>b</sub>**  
MEDICAL DEVICE

07/04/2026 02:16:47 – For information purposes only, can be subject to change

# EAR REUSABLE CLIP SPO2 SENSOR FOR VITALIMATE MINI

Ref. [MONI010E](#)

The ear clip for Mtalimate MINI of Spengler

## Description

Ear clip SpO2 sensor

To be sued with Vitalimate MINI monitor of Spengler.

- length 1m
- 6 pin sensor

Also available for Vitalimate MINI, are the silicon and rigid reusable sensors in every sizes: adult, peditaric and neonat.

## Features

CE Marking	CE0123	Medical device	Class IIb
Bar code	6948309002568	HS Code	9018193090
Warranty (months)	6	Country	China

### DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.055
Bare Unit / Weight (kg)	0.05		

### SPECIFIC STORAGE CONDITIONS

Minimum temperature	-20	Max temperature	55
Minimum humidity	10	Maximum humidity	93

### OTHER FEATURES

📅 07/04/2026 02:16:47 – For information purposes only, can be subject to change