

Ref. MONI010



CE  
0123

DISPOSITIF MEDICAL  
CLASSE  
**II<sub>b</sub>**  
MEDICAL DEVICE

## Description

07/04/2026 01:09:24 – For information purposes only, can be subject to change

=====

**Ref. MONI010**

## Features

CE Marking	CE0123	Medical device	Class IIb
Bar code	6948309002322	HS Code	9018193090
Warranty (months)	6	Country	China

### DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.055
Bare Unit / Weight (kg)	0.05		

### SPECIFIC STORAGE CONDITIONS

Minimum temperature	-20	Max temperature	55
Minimum humidity	10	Maximum humidity	93

### OTHER FEATURES

 07/04/2026 01:09:24 – For information purposes only, can be subject to change