

Ref. BM102S



Description

09/04/2026 23:52:42 – For information purposes only, can be subject to change

=====

Ref. BM102S

Features

CE Marking	CE	Medical device	Class I
Bar code	3700468215738	HS Code	9018901000
Warranty (months)	12		

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	nc

SPECIFIC STORAGE CONDITIONS

Minimum temperature	10	Max temperature	70
Minimum humidity	10	Maximum humidity	100

OTHER FEATURES

 09/04/2026 23:52:42 – For information purposes only, can be subject to change