

**Ref. BM103S**



## Description

09/04/2026 23:37:09 – For information purposes only, can be subject to change

=====

**Ref. BM103S**

## Features

CE Marking	CE	Medical device	Class I
Bar code	3700468215745	HS Code	9018901000
Warranty (months)	12		
<b>DIMENSIONS AND WEIGHT</b>			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	nc
<b>SPECIFIC STORAGE CONDITIONS</b>			
Minimum temperature	10	Max temperature	70
Minimum humidity	10	Maximum humidity	100

### OTHER FEATURES

 09/04/2026 23:37:09 – For information purposes only, can be subject to change