


Ref. TI50004S



CE
0459

GARANTIE
3
3 ANS
GUARANTEE



09/04/2026 23:17:34 – For information purposes only, can be subject to change

=====

Ref. TI50004S



=====

Description

=====

Features

CE Marking	CE0459	Medical device	No
Bar code	3700468215523	HS Code	90189010
Warranty (months)	36	Country	France
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	7	Packed unit / Height (cm)	20
Packed unit / Length (cm)	13	Packed unit / Weight	nc
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	10	Max temperature	70
Minimum humidity	10	Maximum humidity	100

OTHER FEATURES

📅 09/04/2026 23:17:34 – For information purposes only, can be subject to change