


SMALL CHILD CUFF (XS) FOR MOBI

Ref. 513711



CE



12/04/2026 19:58:56 – For information purposes only, can be subject to change

SMALL CHILD CUFF (XS) FOR MOBI

Ref. 513711



Description

Sizes available:

- XS arms with a circumference of 9 to 16 cm (513710)
- S arms from 15 to 26 cm in circumference (513711)

Features

CE Marking	CE	Medical device	Class I
Bar code	3700446028008	HS Code	90189010
Warranty (months)	12	Country	China
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.154
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1	US Size (XXS - XXL)	S

📅 12/04/2026 19:58:56 – For information purposes only, can be subject to change