

**F.O. MACINTOSH BLADE WITH INTEGRATED
LIGHT GUIDE, LUXAMED®,**

Ref. E1.410.012

LUXAMED®



CE



📅 12/04/2026 18:55:30 – For information purposes only, can be subject to change

F.O. MACINTOSH BLADE WITH INTEGRATED LIGHT GUIDE, LUXAMED®,

LUXAMED®Ref. E1.410.012

LUXAMED MACINTOSH laryngoscope blades with integrated fibre optics.

Description

Laryngoscope blade with integrated optical glass

- High-quality stainless steel construction for a long service life
- No dust collecting cracks – easy cleaning and sterilisation
- Maintenance-free
- 8.000 fibres to achieve the maximum service life and illumination
- Optimal light emission due to the largest possible cross section (3 x 7 mm)
- Brightness approx. 14.000 Lux at 3.7 V LED and 2.5 V LED
- Compatible with all handles according to ISO 7376 (green standard)
- Autoclavable up to 134°C – approx. 4.000 times
- Data-Matrix with most important information

Differents size available :

0 : E1.410.012

1 : E1.411.012

2 : E1.412.012

3 : E1.413.012

4 : E1.414.012

5 : E1.415.012

 12/04/2026 18:55:30 – For information purposes only, can be subject to change

**F.O. MACINTOSH BLADE WITH INTEGRATED
LIGHT GUIDE, LUXAMED®,**

LUXAMED®

Ref. E1.410.012

Features

CE Marking	CE	Medical device	Class I
Bar code	04260381960492	HS Code	90189084
Warranty (months)	24	Country	Germany
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	nc
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	5	Max temperature	40
Minimum humidity	10	Maximum humidity	90
OTHER FEATURES			
Size	Size 0		

📅 12/04/2026 18:55:30 – For information purposes only, can be subject to change