

**F.O. MACINTOSH BLADE WITH INTEGRATED  
LIGHT GUIDE, LUXAMED®,**

**LUXAMED®**

Ref. E1.414.012



12/04/2026 20:09:31 – For information purposes only, can be subject to change

# F.O. MACINTOSH BLADE WITH INTEGRATED LIGHT GUIDE, LUXAMED<sup>®</sup>,

**LUXAMED<sup>®</sup>**

Ref. E1.414.012

LUXAMED MACINTOSH laryngoscope blades with integrated fibre optics.

## Description

Laryngoscope blade with integrated optical glass

- High-quality stainless steel construction for a long service life
- No dust collecting cracks – easy cleaning and sterilisation
- Maintenance-free
- 8.000 fibres to achieve the maximum service life and illumination
- Optimal light emission due to the largest possible cross section (3 x 7 mm)
- Brightness approx. 14.000 Lux at 3.7 V LED and 2.5 V LED
- Compatible with all handles according to ISO 7376 (green standard)
- Autoclavable up to 134°C – approx. 4.000 times
- Data-Matrix with most important information

Differents size available :

0 : E1.410.012

1 : E1.411.012

2 : E1.412.012

3 : E1.413.012

4 : E1.414.012

5 : E1.415.012

 12/04/2026 20:09:31 – For information purposes only, can be subject to change

**F.O. MACINTOSH BLADE WITH INTEGRATED  
LIGHT GUIDE, LUXAMED®,**

**LUXAMED®**

Ref. E1.414.012

**Features**

CE Marking	CE	Medical device	Class I
Bar code	04260381960546	HS Code	90189084
Warranty (months)	24	Country	Germany
<b>DIMENSIONS AND WEIGHT</b>			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	nc
<b>SPECIFIC STORAGE CONDITIONS</b>			
Minimum temperature	5	Max temperature	40
Minimum humidity	10	Maximum humidity	90
<b>OTHER FEATURES</b>			
Size	Size 4		

📅 12/04/2026 20:09:31 – For information purposes only, can be subject to change