

# LERICHE RETRACTOR - 28 CM - WIDTH 45 MM

Ref. IV05525



 12/04/2026 20:33:09 – For information purposes only, can be subject to change

# LERICHE RETRACTOR - 28 CM - WIDTH 45 MM



Ref. IV05525

## Features

CE Marking	CE	Medical device	Class I
Bar code	3661040088150	HS Code	90189084
Warranty (months)	24	Country	Pakistan
<b>DIMENSIONS AND WEIGHT</b>			
Packed unit / Width (cm)	21	Bare Unit / Width (cm)	21.6
Packed unit / Height (cm)	4	Bare Unit / Height (cm)	4.6
Packed unit / Length (cm)	28	Bare Unit / Length (cm)	28.1
Packed unit / Weight	0.84	Bare Unit / Weight (kg)	0.401
<b>SPECIFIC STORAGE CONDITIONS</b>			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
<b>OTHER FEATURES</b>			
Packed unit / Quantity	1		

📅 12/04/2026 20:33:09 – For information purposes only, can be subject to change