

EMERGENCY EMERGENCY BAG LIGHT (22 L / 44X25X27CM) - BLUE (EM13.014)

EMERGENCY'S™

Ref. MEEMER02_N



GARANTIE
2
ANS
2 YEARS
GUARANTEE



12/04/2026 20:11:10 – For information purposes only, can be subject to change

EMERGENCY EMERGENCY BAG LIGHT (22 L / 44X25X27CM) - BLUE (EM13.014)

EMERGENCY'S™

Ref. [MEEMER02_N](#)

Practical and lightweight first aid bag with two adjustable internal dividers and a zippered pocket in the flap to facilitate the organization of sanitary equipment

Description

- It also has three outer pockets; the one at the front is divided with rubberized elastic bands
- High resistance zip closures
-
- Reflective strips- Washable case
- Reinforced handles and shoulder straps
- Dimensions 44 x 25 x 27 cm, capacity 22 L
- Material 600D Polyester
- Weight 1.29 kg
-
- Available in waterproof version
-
- Sold empty
- Available in blue (MEEMER02_B), yellow (MEEMER02_J), black (MEEMER02_N), orange (MEEMER02_O), pink (MEEMER02_RC) and red (MEEMER02_R).

 12/04/2026 20:11:10 – For information purposes only, can be subject to change

EMERGENCY EMERGENCY BAG LIGHT (22 L / 44X25X27CM) - BLUE (EM13.014)



Ref. [MEEMER02_N](#)

Features

CE Marking	NONE	Medical device	No
Bar code	8436564542977	HS Code	420299000
Warranty (months)	24	Country	China
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	24	Bare Unit / Width (cm)	26
Packed unit / Height (cm)	22	Bare Unit / Height (cm)	19
Packed unit / Length (cm)	30	Bare Unit / Length (cm)	31.6
Packed unit / Weight	1.3	Bare Unit / Weight (kg)	1.29
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Colors	Black	Packed unit / Quantity	1

📅 12/04/2026 20:11:10 – For information purposes only, can be subject to change