

EMERGENCY EMERGENCY BAG LIGHT (22 L / 44X25X27CM) - BLUE (EM13.014)

EMERGENCY'S™

Ref. MEEMER02_RC



GARANTIE
2 ANS
2
2 YEARS
GUARANTEE



📅 12/04/2026 20:06:35 – For information purposes only, can be subject to change

EMERGENCY EMERGENCY BAG LIGHT (22 L / 44X25X27CM) - BLUE (EM13.014)

EMERGENCY'S™

Ref. [MEEMER02_RC](#)

Practical and lightweight first aid bag with two adjustable internal dividers and a zippered pocket in the flap to facilitate the organization of sanitary equipment

Description

- It also has three outer pockets; the one at the front is divided with rubberized elastic bands
- High resistance zip closures

-

- Reflective strips- Washable case
- Reinforced handles and shoulder straps
- Dimensions 44 x 25 x 27 cm, capacity 22 L
- Material 600D Polyester
- Weight 1.29 kg

-

Available in waterproof version

- Sold empty

- Available in blue (MEEMER02_B), yellow (MEEMER02_J), black (MEEMER02_N), orange (MEEMER02_O), pink (MEEMER02_RC) and red (MEEMER02_R).

EMERGENCY EMERGENCY BAG LIGHT (22 L / 44X25X27CM) - BLUE (EM13.014)

EMERGENCY'S™

Ref. MEEMER02_RC

Features

CE Marking	NONE	Medical device	No
Bar code	8436564542670	HS Code	420299000
Warranty (months)	24	Country	China
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	24	Bare Unit / Width (cm)	25
Packed unit / Height (cm)	22	Bare Unit / Height (cm)	27
Packed unit / Length (cm)	30	Bare Unit / Length (cm)	44
Packed unit / Weight	1.3	Bare Unit / Weight (kg)	1.29
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Colors	Pink	Packed unit / Quantity	1

📅 12/04/2026 20:06:35 – For information purposes only, can be subject to change