

LUXASCOPE OPHTHALMOSCOPE LED 2.5 V - LUXAMED

Ref. F1.416.114

LUXAMED®



CE



12/04/2026 20:31:29 – For information purposes only, can be subject to change

LUXASCOPE OPHTHALMOSCOPE LED 2.5 V - LUXAMED

LUXAMED[®]

Ref. F1.416.114

Direct ophthalmoscope in pocket size

Description

Direct ophthalmoscope in pocket size

Refraction balance: + 20 bis - 20 D (\pm 1 2 3 4 6 8 10 15 20 D)

5 different apertures: reticle, large circle, small circle, semicircle, green filter

Illumination: LED

Can be dimmed for an individual glare-free brightness

Automatic switchoff after 3 minutes

Life expectancy of LEDs: approx. 50.000 hours by constant current drive

Colour temperature: approx. 4.200 K

Optimized design of head

Two-component handle (aluminium/fibre-glass reinforced), also combinable with LuxaScope LED otoscope- and dermatoscope head

Electronic ON/OFF push-button switch

Power supply: 2x AAA Alkaline batteries, inclusive

Available in :

- White : C1.516.914
- Blue : C1.516.934
- Pink : C1.516.974
- Yellow : C1.516.944
- Kids : C1.916.914
- Black : C1.516.114
- Green : C1.516.954
- Purple : C1.516.964

 12/04/2026 20:31:29 – For information purposes only, can be subject to change

**LUXASCOPE OPHTHALMOSCOPE LED 2.5 V
- LUXAMED**

LUXAMED®

Ref. F1.416.114

Features

CE Marking	CE	Medical device	Class I
Bar code	04260381964315	HS Code	90185090
Warranty (months)	24	Country	Germany
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	nc
Bare Unit / Weight (kg)	0.12		
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	5	Max temperature	40
Minimum humidity	10	Maximum humidity	90
OTHER FEATURES			
Colors	Black		

 12/04/2026 20:31:29 – For information purposes only, can be subject to change