

Ref. 335111



Description

Cuffs with dual tube bladder & coupling (60 cm tube with Classic® bulb and valve, and 10 cm with male screw-on metal fitting)

Available in :

Adult (M), nylon, velcro, black (335111)

Adult (L), nylon, black (335101)

15/04/2026 00:52:54 – For information purposes only, can be subject to change

Ref. 335111



Features

CE Marking	CE	Medical device	Class I
Bar code	3700446028435	HS Code	9018901000
Warranty (months)	6	Country	Indonesia
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.245
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

15/04/2026 00:52:54 – For information purposes only, can be subject to change