

MAGISTER STETHOSCOPE FOR ADULT / BLACK

Ref. 507200



16/04/2026 00:18:12 – For information purposes only, can be subject to change

MAGISTER STETHOSCOPE FOR ADULT / BLACK

Ref. 507200



Description

For general-purpose auscultation. Stainless steel dual head chestpiece for adults. Multifrequency diaphragm. Angled binaural with built-in spring and flexible, screw-on tips. Ideal grip and high performance make it the perfect product for general and specialist use.

Available in en Grey (ref 507204), Red (ref: 507201), Green (ref :507209), Bleu (ref : 507202), Black (ref: 507200)

STORAGE CONDITIONS :

Atmospheric pressure : 700-1060 hPa

Temperature: -10° to 40° C

Humidity: 95% max

 16/04/2026 00:18:12 – For information purposes only, can be subject to change

MAGISTER STETHOSCOPE FOR ADULT / BLACK



Ref. 507200

Features

CE Marking	CE	Medical device	No
Bar code	3700446009069	HS Code	90181990
Warranty (months)	60		
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	16	Packed unit / Height (cm)	4
Packed unit / Length (cm)	32	Packed unit / Weight	0.398
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	-10	Max temperature	40
Minimum humidity	95	Maximum humidity	95
OTHER FEATURES			
Colors	Black	Packed unit / Quantity	1

16/04/2026 00:18:12 – For information purposes only, can be subject to change