


# MAGISTER STETHOSCOPE FOR ADULT / BLACK

Ref. 507209



16/04/2026 00:43:04 – For information purposes only, can be subject to change

# MAGISTER STETHOSCOPE FOR ADULT / BLACK

Ref. 507209



## Description

For general-purpose auscultation. Stainless steel dual head chestpiece for adults. Multifrequency diaphragm. Angled binaural with built-in spring and flexible, screw-on tips. Ideal grip and high performance make it the perfect product for general and specialist use.

Available in en Grey ( ref 507204), Red (ref: 507201), Green (ref :507209), Bleu (ref : 507202), Black (ref: 507200)

### STORAGE CONDITIONS :

Atmospheric pressure : 700-1060 hPa

Temperature: -10° to 40° C

Humidity: 95% max

# MAGISTER STETHOSCOPE FOR ADULT / BLACK



Ref. 507209

## Features

|                                    |               |                           |                 |
|------------------------------------|---------------|---------------------------|-----------------|
| CE Marking                         | CE            | Medical device            | No              |
| Bar code                           | 3700446004002 | HS Code                   | 9018.90.20.00-9 |
| Warranty (months)                  | 60            |                           |                 |
| <b>DIMENSIONS AND WEIGHT</b>       |               |                           |                 |
| Packed unit / Width (cm)           | nc            | Packed unit / Height (cm) | nc              |
| Packed unit / Length (cm)          | nc            | Packed unit / Weight      | 0.398           |
| <b>SPECIFIC STORAGE CONDITIONS</b> |               |                           |                 |
| Minimum temperature                | -10           | Max temperature           | 40              |
| Minimum humidity                   | 95            | Maximum humidity          | 95              |
| <b>OTHER FEATURES</b>              |               |                           |                 |
| Colors                             | Green         | Packed unit / Quantity    | 1               |

16/04/2026 00:43:04 – For information purposes only, can be subject to change