


STETHOSCOPE MAGISTER II DUAL, SILVER

Ref. 507252



CE



16/04/2026 00:42:15 – For information purposes only, can be subject to change

STETHOSCOPE MAGISTER II DUAL, SILVER

Ref. 507252



The Magister II is a high-performance professional stethoscope for versatile auscultation.

Description

HAUTE QUALITÉ ACOUSTIQUE

- Membrane flottante multifréquences
- Pression légère = écoute des basses fréquences (vasculaires et cardiaques)
- Pression ferme = écoute des hautes fréquences (pulmonaires)
- Tubulure simple avec un conduit haute performance
- Lyres en aluminium

CONFORT

- Bague antifroid
- Embouts auriculaires en silicone
- Nouveau design de lyres plus adapté à l'anatomie
- Indicateur de pavillon

Disponible en version Silver (507251), Black (507261), Myrtille (507262) et Corail (507252)

CONDITIONS DE STOCKAGE :

Pression atmosphérique : 700-1060 hPa

Température : -10° à 40° C

Humidité : 95% max

STETHOSCOPE MAGISTER II DUAL, SILVER

Ref. 507252



Features

CE Marking	CE	Medical device	Class I
Bar code	3700446042134	HS Code	90181990
Warranty (months)	60		
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	16	Packed unit / Height (cm)	4
Packed unit / Length (cm)	32	Packed unit / Weight	0.43
Bare Unit / Weight (kg)	0.105		
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	-10	Max temperature	40
Minimum humidity	NA	Maximum humidity	95
OTHER FEATURES			
Colors	Coral	Packed unit / Quantity	1

📅 16/04/2026 00:42:15 – For information purposes only, can be subject to change