

IUD (STERILET) IMPLANTATION AND REMOVING SURGICAL SET

Ref. IBOI118



15/04/2026 00:07:07 – For information purposes only, can be subject to change

IUD (STERILET) IMPLANTATION AND REMOVING SURGICAL SET



Ref. IBOI118

Description

This box contains:

- Stainless steel box, 35 × 18 × 8 cm – IS
- Spouted cup, stainless steel, Ø 80 mm, 165 ml – IS
- Rochester-Péan / Kocher forceps, straight/jaws serrated, curved, 22 cm
- Pozzi forceps, straight, 24 cm – IS, Holtex
- Foerster forceps, serrated, straight, 24 cm, Holtex
- Novak biopsy cannula, 23 cm (curette probe)
- Circéron forceps for IUD grasping, 28 cm – IS
- Sims scissors, blunt/blunt, curved, 20 cm – IS, Holtex
- Valleix uterine sound with slider, 33 cm
- Cusco vaginal speculum, folding, 32 mm, length 90 mm
- Dressing forceps, long pattern, straight, 24 cm – IS, Holtex
- Cotton holder (applicator), 30 cm
- IUD removal hook, 25 cm
- Instrument holder, Ø 52 mm, height 200 mm
- Therrum forceps, 27 cm

Features

CE Marking	CE	Medical device	Class I
Bar code	3661040081755	HS Code	90189084
Warranty (months)	24	Country	Pakistan
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	2.629
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

📅 15/04/2026 00:07:07 – For information purposes only, can be subject to change