


**Ref. TS53000S**



**CE**  
**0459**



📅 15/04/2026 23:57:27 – For information purposes only, can be subject to change

=====

**Ref. TS53000S**



=====

## Description

=====

## Features

CE Marking	CE0459	Medical device	No
Bar code	3700468215486	HS Code	90189010
Country	France		
<b>DIMENSIONS AND WEIGHT</b>			
Packed unit / Width (cm)	7	Packed unit / Height (cm)	25
Packed unit / Length (cm)	16	Packed unit / Weight	nc
<b>SPECIFIC STORAGE CONDITIONS</b>			
Minimum temperature	10	Max temperature	70
Minimum humidity	10	Maximum humidity	100

### OTHER FEATURES

📅 15/04/2026 23:57:27 – For information purposes only, can be subject to change