

SUNTECH MONOBLOC DOUBLE TUBE CUFF REUSABLE EXTRA LARGE (38-50CM) - WITHOUT FITTING

SUNTECH

Ref. DSU231049



CE



Description

Suntech brown nylon armband individually.
Double tubing, without fitting, wide size adult 38-50 cm

16/04/2026 23:45:23 – For information purposes only, can be subject to change

**SUNTECH MONOBLOC DOUBLE TUBE CUFF
REUSABLE EXTRA LARGE (38-50CM) -
WITHOUT FITTING**

SUNTECH

Ref. DSU231049

Features

CE Marking	CE	Medical device	Class I
Bar code	3700446013998	HS Code	9018905040
Warranty (months)	24	Country	China
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.769
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

 16/04/2026 23:45:23 – For information purposes only, can be subject to change