

**SET OF 5 SUNTECH MONOBLOC REUSABLE
ARMBANDS EXTRA LARGE (38-50CM) -
BAYONET METAL CONNECTOR**

SUNTECH

Ref. DSU231116



CE



Intended for non-invasive blood pressure measurement in adults.

Description

Reusable nylon cuff sold in boxes of 5 units.
Single tubing, bayonette connection 38-50cm.

16/04/2026 23:43:28 – For information purposes only, can be subject to change

SET OF 5 SUNTECH MONOBLOC REUSABLE ARMBANDS EXTRA LARGE (38-50CM) - BAYONET METAL CONNECTOR

SUNTECH

Ref. DSU231116

Features

CE Marking	CE	Medical device	Class I
Bar code	3661041008843	HS Code	9018905040
Warranty (months)	24	Country	China
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	25	Packed unit / Height (cm)	10
Packed unit / Length (cm)	23	Packed unit / Weight	0.8
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

16/04/2026 23:43:28 – For information purposes only, can be subject to change