

## COOL BAG FOR DIABETIQUE PM DIA, SILVER

Ref. MEDIA02\_A

**ELITE BAGS**



16/04/2026 23:42:51 – For information purposes only, can be subject to change

# COOL BAG FOR DIABETIQUE PM DIA, SILVER



Ref. [MEDIA02\\_A](#)

**Kit designed with the spectrum of patients with diabetes in mind**

## Description

- Double insulated compartment to keep products cool
- Includes reusable cold gel
- Material: 420D fabric
- Weight: 130 g
- Exterior dimensions: 19 x 10 x 5 cm / 0.95 l
- Max. recommended capacity: 600 g
- Available in blue (ref: MEDIA02\_B) or silver (ref: MEDIA02\_A)
- Sold empty

## Features

|                   |               |                |               |
|-------------------|---------------|----------------|---------------|
| CE Marking        | NONE          | Medical device | No            |
| Bar code          | 8436564540744 | HS Code        | 4202.12.99.90 |
| Warranty (months) | 24            | Country        | China         |

### DIMENSIONS AND WEIGHT

|                           |      |                         |      |
|---------------------------|------|-------------------------|------|
| Packed unit / Width (cm)  | 20   | Bare Unit / Width (cm)  | 10   |
| Packed unit / Height (cm) | 3    | Bare Unit / Height (cm) | 5    |
| Packed unit / Length (cm) | 14   | Bare Unit / Length (cm) | 19   |
| Packed unit / Weight      | 0.25 | Bare Unit / Weight (kg) | 0.13 |

### SPECIFIC STORAGE CONDITIONS

|                     |    |                  |    |
|---------------------|----|------------------|----|
| Minimum temperature | NA | Max temperature  | NA |
| Minimum humidity    | NA | Maximum humidity | NA |

### OTHER FEATURES

|        |         |                        |   |
|--------|---------|------------------------|---|
| Colors | Argenté | Packed unit / Quantity | 1 |
|--------|---------|------------------------|---|

📅 16/04/2026 23:42:51 – For information purposes only, can be subject to change