

**Ref. 334129**



CE



19/04/2026 21:10:00 – For information purposes only, can be subject to change

Ref. **334129**


## Description

Nylon cuffs with quick coupling for Lian® with 60 cm single tubing bladder fitted with a quick coupling for Lian® NM, Lian® Métal, Lian® Nano, Lian® Scénic, Ergo Sphyg®, Lian® Color. Washable at 30-40° C.

## Features

|                                    |               |                           |            |
|------------------------------------|---------------|---------------------------|------------|
| CE Marking                         | CE            | Medical device            | Class I    |
| Bar code                           | 3700446000011 | HS Code                   | 9018901000 |
| Warranty (months)                  | 6             | Country                   | Indonesia  |
| <b>DIMENSIONS AND WEIGHT</b>       |               |                           |            |
| Packed unit / Width (cm)           | nc            | Packed unit / Height (cm) | nc         |
| Packed unit / Length (cm)          | nc            | Packed unit / Weight      | 0.11       |
| <b>SPECIFIC STORAGE CONDITIONS</b> |               |                           |            |
| Minimum temperature                | NA            | Max temperature           | NA         |
| Minimum humidity                   | NA            | Maximum humidity          | NA         |
| <b>OTHER FEATURES</b>              |               |                           |            |
| Packed unit / Quantity             | 1             |                           |            |

📅 19/04/2026 21:10:00 – For information purposes only, can be subject to change