

CHARRIER SAW 27 CM

Ref. IS15227



CE



- Charrier amputation saw, 27 cm long, designed for surgical procedures.
- Durable serrated blade with an ergonomic handle ensuring a secure and controlled grip.

📅 20/04/2026 22:55:25 – For information purposes only, can be subject to change

CHARRIER SAW 27 CM

Ref. IS15227



Features

CE Marking	CE	Medical device	Class Ir
Bar code	3661040030845	HS Code	90189084
Warranty (months)	24	Country	Pakistan
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	10	Packed unit / Height (cm)	5
Packed unit / Length (cm)	31	Bare Unit / Length (cm)	27
Packed unit / Weight	0.247	Bare Unit / Weight (kg)	0.247
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

📅 20/04/2026 22:55:25 – For information purposes only, can be subject to change