

SOFT CUFF® NIBP CUFF – 5 SIZES

Ref. BT105S

CE



22/04/2026 22:25:13 – For information purposes only, can be subject to change

SOFT CUFF® NIBP CUFF – 5 SIZES

Ref. BT105S

The ultra-comfortable Soft Cuff® cuffs!

Description

The Soft Cuff® cuff with pocket offers comfort and softness for the patient.

Suitable for all types of patients, with 6 sizes available (5 cuffs and one thigh cuff)

Made from Oxford nylon, a material that feels pleasant against the skin

Latex-free

Compatible with the Lian Dessillons SPENGLER, COMEN monitors and others.

Available in 5 sizes:

BT101S infant size

BT102S child size

BT103S small adult size

BT104S adult size

BT105S large adult size

SPENGLER MEDICAL IS COMMITTED TO SUSTAINABILITY!

For Soft Cuff® arm and leg cuffs, SPENGLER MEDICAL has replaced plastic packaging with FSC-certified kraft envelopes.

SOFT CUFF® NIBP CUFF – 5 SIZES

Ref. BT105S

Features

CE Marking	CE	Medical device	Class I
Bar code	3700468215707	HS Code	9018901000
Warranty (months)	12	Country	France

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	nc

SPECIFIC STORAGE CONDITIONS

Minimum temperature	10	Max temperature	70
Minimum humidity	10	Maximum humidity	100

OTHER FEATURES

 22/04/2026 22:25:13 – For information purposes only, can be subject to change