

## REUSABLE NIBP FOR VITALIMATE MINI

Ref. MONI011G



CE



📅 22/04/2026 02:07:11 – For information purposes only, can be subject to change

## REUSABLE NIBP FOR VITALIMATE MINI

Ref. [MONI011G](#)

Vitalimate MINI reusable NIBP cuffs

### Description

The Vitalimate MINI NIBP cuffs are comfortable and come in all sizes for all types of patients.

- Reusable adult NIBP cuff (27-35cm) : ref MONI011
- Reusable large adult NIBP cuff (34-43cm) : ref MONI011A
- Small reusable adult NIBP cuff (20.5-28 cm): ref MONI011G
- Reusable pediatric NIBP cuff (14-21.5 cm): ref MONI011C
- Reusable neonatal NIBP cuff (10-15 cm): ref MONI011D
- Reusable adult NIBP thigh cuff (42-54 cm): ref MONI011B

Disposable neonatal cuffs are also available for infants.

### Features

CE Marking	CE	Medical device	Class I
Bar code	6948309002575	HS Code	9018193090
Warranty (months)	6	Country	China

#### DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.085
Bare Unit / Weight (kg)	0.08		

#### SPECIFIC STORAGE CONDITIONS

Minimum temperature	-20	Max temperature	55
Minimum humidity	10	Maximum humidity	93

#### OTHER FEATURES

22/04/2026 02:07:11 – For information purposes only, can be subject to change