


**POLYESTER AUTOCLAVABLE CUFF
STRAIGHT WITH 1 BLADDER FOR DDM
PNEUMATIC TOURNIQUET**

Ref. GBS102

DDM 



CE



📅 23/04/2026 23:05:03 – For information purposes only, can be subject to change

POLYESTER AUTOCLAVABLE CUFF STRAIGHT WITH 1 BLADDER FOR DDM PNEUMATIC TOURNIQUET



Ref. GBS102

Those cuffs are autoclavable!

Description

Cuffs for DDM pneumatic tourniquets are made of polyester. They are disinfectable and autoclavable and composed of a removable silicon bladder for easy autoclaving.

Those ones are designed in straight form

Several sizes are available

Dimensions of the bladder in cm

Cuff infant (GBS101)	19 x 4 cm
Cuff child (GBS102)	26.5 x 5 cm
Cuff small adult (GBS101)	35.5 x 6 cm
Cuff adult (GBS103)	53 x 6 cm

**POLYESTER AUTOCLAVABLE CUFF
STRAIGHT WITH 1 BLADDER FOR DDM
PNEUMATIC TOURNIQUET**



Ref. GBS102

Features

CE Marking	CE	Medical device	Class I
Bar code	3700468203483	HS Code	90189084
Warranty (months)	24	Country	France
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	nc
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	-10	Max temperature	40
Minimum humidity	30	Maximum humidity	40

OTHER FEATURES

23/04/2026 23:05:03 – For information purposes only, can be subject to change