

REUSABLE SPO2 SENSOR ADU FINGER-CLIP (RIGID) FOR NC6 NC

COMEN

Ref. [BCCAPT09](#)

The reusable rigid sensor for Comen monitors

Description

The reusable rigid clip adult SpO2 sensor for Comen monitors.

Compatible with NC6, NC7 and ND10

Features

CE Marking	CE1639	Medical device	Class IIb
Bar code	6945429080636	HS Code	9018193090
Country	China		

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	14	Packed unit / Height (cm)	1
Packed unit / Length (cm)	21	Packed unit / Weight	0.065
Bare Unit / Weight (kg)	0.045		

SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

OTHER FEATURES

26/04/2026 19:24:23 – For information purposes only, can be subject to change