

**REUSABLE SPO2 WRAP SENSOR
PEDIATRIC AND NEONAT FOR COMEN
MONITORS**

COMEN

Ref. BCCAPT11



CE
1639

DISPOSITIF MEDICAL
CLASSE
II_b
MEDICAL DEVICE

📅 26/04/2026 18:36:49 – For information purposes only, can be subject to change

REUSABLE SPO2 WRAP SENSOR PEDIATRIC AND NEONAT FOR COMEN MONITORS

COMEN

Ref. BCCAPT11

The SpO2 sensor for the little ones!

Description

Reusable SpO2 wrap sensor for pediatric and neonat patient to be used with Comen monitors.

Compatible with NC6, NC7 and ND10

Features

CE Marking	CE1639	Medical device	Class IIb
Bar code	6945429080742	HS Code	9018193090
Country	China		

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	nc

SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

OTHER FEATURES

📅 26/04/2026 18:36:49 – For information purposes only, can be subject to change