

Ref. 336401



CE



Description

Nylon cuffs with quick coupling for Lian® with 60 cm single tubing bladder fitted with a quick coupling for Lian® NM, Lian® Métal, Lian® Nano, Lian® Scénic, Ergo Sphyg®, Lian® Color. Washable at 30-40° C.

27/04/2026 23:09:26 – For information purposes only, can be subject to change

=====

Ref. 336401



Features

CE Marking	CE	Medical device	Class I
Bar code	3700446026615	HS Code	9018901000
Warranty (months)	6	Country	China
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.29
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

 27/04/2026 23:09:26 – For information purposes only, can be subject to change