

IDEAL+ SINGLE HEAD ADULT, BLUE

Ref. DH7081_N



27/04/2026 23:04:58 – For information purposes only, can be subject to change

IDEAL+ SINGLE HEAD ADULT, BLUE

Ref. DH7081_N



Single-head adult stethoscope, designed for reliable and accurate auscultation.

Description

- Stethoscope with single adult metal horn (Ø: 45 mm)•
- Anti-cold ring and horn of the color of the tube•
- Aluminum lyre incorporated into the tubing
- Also available in double pavilion for adults and a pediatric model

Available in Black (DH7081_N) or Red (DH7081_R)

Features

CE Marking	CE	Medical device	Class I
Bar code	3661041000120	HS Code	9018908400
Warranty (months)	12	Country	China

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.15

SPECIFIC STORAGE CONDITIONS

Minimum temperature	-20	Max temperature	70
Minimum humidity	20	Maximum humidity	85

OTHER FEATURES

Colors	Black	Packed unit / Quantity	1
--------	-------	------------------------	---

📅 27/04/2026 23:04:58 – For information purposes only, can be subject to change