

**LIAN METAL / GREY NYLON CUFF / OBESE  
(XL)**

**Ref. 521330**



**CE**  
**0459**

DISPOSITIF MEDICAL  
CLASSE  
**Im**  
MEDICAL DEVICE


GARANTIE  
**3**  
3 ANS  
3 YEARS  
GUARANTEE

29/04/2026 00:03:01 – For information purposes only, can be subject to change

# LIAN METAL / GREY NYLON CUFF / OBESE (XL)



**Ref. 521330**



## Description

Large dial, ambidextrous model (simply rotate to change hands). Metal casing with large dial for easy reading from all angles. Soft PVC manometer cover. Ergonomic spoon handle and inflation bulb. Cuffs with screw-on metal quick coupling. Supplied in a nylon pouch.

Available in :

- Adult (M), cotton, velcro, red (521317)
- Adult (M), cotton, velcro, grey (521311)
- Adult large (L), cotton, velcro, grey (521310)
- Child/Small Adult (S), nylon, velcro, grey (521316)
- Multibrassard (S, M, L) velcro (521330)

## Features

CE Marking	CE0459	Medical device	Class Im (measuring function)
Bar code	3700446009939	HS Code	90189010
Warranty (months)	36	Country	France
<b>DIMENSIONS AND WEIGHT</b>			
Packed unit / Width (cm)	18	Packed unit / Height (cm)	9
Packed unit / Length (cm)	26	Packed unit / Weight	1.213
<b>SPECIFIC STORAGE CONDITIONS</b>			
Minimum temperature	-20	Max temperature	70
Minimum humidity	85	Maximum humidity	85
<b>OTHER FEATURES</b>			
Packed unit / Quantity	1		

📅 29/04/2026 00:03:01 – For information purposes only, can be subject to change