


# VAQUEZ-LAUBRY CLASSIC VELCRO CUFF WITH BLADDER/DARK BLUE/LARGE ADULT



Ref. 330011



30/04/2026 01:13:43 – For information purposes only, can be subject to change

# VAQUEZ-LAUBRY CLASSIC VELCRO CUFF WITH BLADDER/DARK BLUE/LARGE ADULT



Ref. 330011

## Description

Made of polycotton with inner bladder and dual tubing (10 cm & 60 cm) - also suited for Vaquez-Laubry® with full-grain leather manometer holder, and Velcro fastening.

Adult (M) cotton, velcro, grey: 330011

Adult (M) cotton, velcro, red : 330013

Adult (M) navy, with straps: 330112

Adult (M) red, with straps: 330113

Adult (L) cotton, velcro, navy: 330002

Child/small adult (S), cotton, velcro, grey: 330021

## Features

CE Marking	NONE	Medical device	No
Bar code	3700446028534	HS Code	9018901000
Warranty (months)	6	Country	France

### DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.195

### SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

### OTHER FEATURES

Packed unit / Quantity	1
------------------------	---

📅 30/04/2026 01:13:43 – For information purposes only, can be subject to change